




# PAD Trolley User Manual

*Model: RS010-100-11*



# CONTENT

1	Safety Information .....	2
2	Equipment Symbols .....	3
3	Area of application/ Intended use .....	3
4	Reporting obligation .....	3
5	CE Mark .....	4
6	EU Authorized Representative .....	4
7	Manual Safekeeping .....	4
8	Specification .....	5
9	Packing List .....	5
10	Roll Stand Parts Schematic .....	6
11	Assemble the roll stand .....	7
12	Install PAD .....	9
13	Maintenance and Adjustment .....	10
14	Cleaning and disinfection .....	10

# 1 Safety Information

## 1.1 Dangers

There are no dangers that refer to the product in general. Specific "Danger" statements may be given in the respective sections of this operation manual.<sup>1</sup>

## 1.2 Warnings

---

### WARNING

---

- **Do not dismantle the wall frame by yourself.**

The maintenance of the roll stand must be carried out by maintenance personnel trained and authorized by the company.

- **Handle packaging materials**

You must comply with local regulations or the hospital's waste disposal system. Packaging materials must be placed out of reach of children.

- **MORTAL DANGER BY ELECTRIC SHOCK**

The electric cables of ME equipment and ME systems may only be connected to the power source and to the medical device respectively after finishing installation. Prior to maintenance or dis-assembly all power cables must be disconnected from the power source.

---

## 1.3 Note

---

### Note

---









- **MOUNTING THE MEDICAL ELECTRICAL EQUIPMENT**

- Only when all screws are tightened and securely fixed can the PAD be mounted.
-

## 2 Equipment Symbols

### NOTE

- Some symbols may not appear on all equipment.

	Attention: Consult accompanying documents (this manual).		Serial number
	Read Operator's Manual		CE mark
	Manufacturer		Importer
	Manufacture date		Medical Device

## 3 Area of application/ Intended use

CORITON Medical Equipment Roll Stand are intended to mount and position medical electrical equipment (hereafter referred to as ME equipment) and accessories to those and are used in operating theaters as well as treatment rooms as auxiliary devices for doctors and medical staff.

## 4 Reporting obligation

Reporting obligation

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **5 CE Mark**



## **6 EU Authorized Representative**

**Umedwings Netherlands B.V.**

ADD: Treubstraat 1, 2288 EG, Rijswijk, The Netherlands

Tel: +31 (0) 642758955 Email: [ar@umedwings.eu](mailto:ar@umedwings.eu)

## **7 Manual Safekeeping**

For any questions that may arise at a later time please safe-keep our manual!

## 8 Specification


### Hardware Specifications

Model	RS010-100-11
Load Capacity	15kg
Lifting method	Mechanical lifting
Trolley Size	600*600*1500mm
Adjustable Height	1250~1500mm
Basket	330*210*130mm, Environment-friendly iron, Stoving varnish-Medical White
Material	Aluminum alloy+Sheet metal
Package Size	78*63*21.5 cm
Net Weight	14kg±0.5kg
Gross Weight	17kg±0.3kg
Certificates	ISO9001/ROSH/FDA/ISO13485




## 9 Packing List

NO.	Description	Qty.	Note
1	Base with caster	1	5 with breaker
2	Roll Stand Pole	1	
3	Basket	1	Contains 2 pairs of lock ring
4	Handle	1	
5	PAD arm	1	
6	#3 Hexagon wrench	1	Install the lock ring and basket screws
7	#4 Hexagon wrench	1	
8	#5 Hexagon wrench	1	Install the metal lock ring
9	#6 Hexagon wrench	1	Install the screws for fixing the base of the vertical pole



## 10 Roll Stand Parts Schematic







# 11 Assemble the roll stand

<div>1</div>	<ul style="list-style-type: none"><li>- Insert pole to base hole.</li><li>- Use 4 PCS screws (M8*40) to fix the pole to the base</li></ul> <p>Note: When the pole is inserted into the base, pay attention to the limit opening between the pole and the base.</p>	
<div>2</div>	<ul style="list-style-type: none"><li>- Use 4 socket head cap screws (M4*14) to fix 2 pairs of lock ring on the pole</li><li>- Use 4 socket head cap screws(M4*10) to fasten the storage basket to the lock ring</li></ul>	
<div>3</div>	<ul style="list-style-type: none"><li>-Use 2 socket head cap screws (M4*10) to fix the handle to the pole</li></ul>	



<p>4</p>	<p>-Use 2 socket head cap screws (M6*30) to fix the support arm on the cart lift rod</p>	
<p>5</p>	<p>-Use 2 socket head cap screws (M5*45 ) to fix the PAD clamp to the support arm</p>	

## 12 Install PAD

1	Installation: adjust the clamp to a fixed position of the pad	
2	Place the plate diagonally on the fixture	
3	Place it and lock the rear key.  Note: Adjust the clamp according to the size of the tablet. As shown in the figure, the screw less end can be pushed or pulled, or the large clamp can be replaced.	
4	Fix the tablet computer on the fixture, pay attention to the diagonal installation, lock the key	

## 13 Maintenance and Adjustment.

Very little maintenance is required during the life span of the roll stand, however it is recommended to check all pivot /swivel points as well as all bearing points regularly . Please check all screws on a regular basis (at least once a year ) and tighten them if necessary.

## 14 Cleaning and disinfection

The roll stand may be cleaned with mild , non-abrasive solutions commonly used in hospital environment (e.go diluted bleach , ammonia , or alcohol solutions).

The surface finish may be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.

Do not use steel wool or other abrasive materials. Avoid all liquids entering the roll stand.

**CAUTION:** Our company makes no claims concerning the efficacy of the disinfection process. Please consult your hospital's infection control manage.